

## Product Datasheet

### Human Hepcidin-25 ELISA Kit (orb407172)

|                          |                                       |
|--------------------------|---------------------------------------|
| <b>Description</b>       | Human Hepcidin<br>25(Hepc25)ELISA Kit |
| <b>Reactivity</b>        | Human                                 |
| <b>Range</b>             | 15.6 ng/ml -1000 ng/ml                |
| <b>Target</b>            | Hepc25                                |
| <b>Note</b>              | For research use only                 |
| <b>Application notes</b> | Detection Wavelength: 450<br>nm       |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates  |
| <b>Sensitivity</b>       | 3.9 ng/ml                             |
| <b>Expiration Date</b>   | Please enquire.                       |