

Product Datasheet

Human CYP2R1 ELISA Kit (orb407209)

Description	Human Vitamin D 25-hydroxylase(CYP2R1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CYP2R1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6VVX0
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.