

## Product Datasheet

### Human MG7-Ag ELISA Kit (orb407213)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Human MG7-Ag ELISA Kit               |
| <b>Reactivity</b>        | Human                                |
| <b>Range</b>             | 1 U/mL-15 U/mL                       |
| <b>Target</b>            | MG7-Ag                               |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450 nm         |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates |
| <b>Sensitivity</b>       | 0.25 U/mL                            |
| <b>Expiration Date</b>   | Please enquire.                      |