

Product Datasheet

Human PYCARD ELISA Kit (orb407233)

Description	Human Apoptosis-associated speck-like protein containing a CARD(PYCARD) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	PYCARD
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9ULZ3
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.