

Product Datasheet

Human KRT18 ELISA Kit (orb407242)

Description	Human cytokeratin 18,CK-18 ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-15 ng/mL
Target	KRT18
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell_x000D_ culture supernates
Uniprot ID	P05783
Sensitivity	0.25 ng/mL
Expiration Date	Please enquire.