

Product Datasheet

Human CYP3A4 ELISA Kit (orb407255)

Description	Human Cytochrome P450 3A4(CYP3A4) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Target	CYP3A4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P08684
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.