

Product Datasheet

Human FNDC5 ELISA Kit (orb407277)

Description	Human Fibronectin type III domain-containing protein 5(FNDC5) ELISA kit
Reactivity	Human
Range	7.8 ng/mL-500 ng/mL
Target	FNDC5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8NAU1
Sensitivity	1.95 ng/mL
Expiration Date	Please enquire.