

Product Datasheet

Human ABCG2 ELISA Kit (orb407281)

Description	Human ATP-binding cassette transporter G2,ABCG2 ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	ABCG2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell lysates, cell culture supernates
Uniprot ID	Q9UNQ0
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.