

## Product Datasheet

### Human TNNT2 ELISA Kit (orb407293)

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Description</b>       | Human cardiac isoform of Troponin T,cTnT ELISA Kit |
| <b>Reactivity</b>        | Human                                              |
| <b>Range</b>             | 15.6 pg/ml-1000 pg/ml                              |
| <b>Target</b>            | TNNT2                                              |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                  |
| <b>Uniprot ID</b>        | <b>P45379</b>                                      |
| <b>Sensitivity</b>       | 3.9 pg/ml                                          |
| <b>Expiration Date</b>   | Please enquire.                                    |