

Product Datasheet

Human IKZF3 ELISA Kit (orb407299)

Description	Human Zinc finger protein Aiolos(IKZF3) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	IKZF3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UKT9
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.