

Product Datasheet

Human AGTR1 ELISA Kit (orb407323)

Description	Human angiotensin II receptor 1,ANGIIR-1 ELISA Kit
Reactivity	Human
Range	46.88 pg/mL-3000 pg/mL
Target	AGTR1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P30556
Sensitivity	11.72 pg/mL
Expiration Date	Please enquire.