

Product Datasheet

Human VEGF121 ELISA Kit (orb407330)

Description	Human Vascular Endothelial Growth Factor 121(VEGF121) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	VEGF121
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P15692
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.