

Product Datasheet

Human FLT1 ELISA Kit (orb407331)

Description	Human vascular endothelial cell growth factor receptor 1 (VEGFR-1/Flt1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	VEGFR-1/Flt1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P17948
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.