

Product Datasheet

Human PDGFC ELISA Kit (orb407335)

Description	Human vascular endothelial cell growth factor E,VEGF-E ELISA Kit
Reactivity	Human
Range	6.25 pg/mL-400 pg/mL
Target	PDGFC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, urine, cell culture supernates, tissue homogenates
Uniprot ID	Q9NRA1
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.