

Product Datasheet

Human PDGFR beta ELISA Kit (orb407385)

Description	Human Beta-type platelet-derived growth factor receptor(PDGFRB) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	PDGFRB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09619
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.