

## Product Datasheet

### Human IGF-1 ELISA Kit (orb407409)

|                          |                                                        |
|--------------------------|--------------------------------------------------------|
| <b>Description</b>       | Human insulin-like growth factors<br>1,IGF-1 ELISA Kit |
| <b>Reactivity</b>        | Human                                                  |
| <b>Range</b>             | 7.8 ng/mL-500 ng/mL                                    |
| <b>Target</b>            | IGF-1                                                  |
| <b>Note</b>              | For research use only                                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                      |
| <b>Uniprot ID</b>        | <b>P05019</b>                                          |
| <b>Sensitivity</b>       | 1.95 ng/mL                                             |
| <b>Expiration Date</b>   | Please enquire.                                        |