

Product Datasheet

Human IGF2 ELISA Kit (orb407412)

Description	Human insulin-like growth factors 2,IGF-2 ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	IGF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, tissue homogenates
Uniprot ID	P01344
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.