

Product Datasheet

Human IGF2BP3 ELISA Kit (orb407413)

Description	Human Insulin-like growth factor 2 mRNA-binding protein 3(IGF2BP3) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Target	IGF2BP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O00425
Sensitivity	9.37 pg/mL
Expiration Date	Please enquire.