

Product Datasheet

Human IGFBP1 ELISA Kit (orb407415)

Description	Human insulin-like growth factors binding protein 1,IGFBP-1 ELISA Kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	IGFBP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P08833
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.