

Product Datasheet

Human IGFBP2 ELISA Kit (orb407416)

Description	Human insulin-like growth factors binding protein 2,IGFBP-2 ELISA Kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	IGFBP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P18065
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.