

Product Datasheet

Human IGFBP5 ELISA Kit (orb407419)

Description	Human insulin-like growth factor binding protein 5,IGFBP-5 ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	IGFBP5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	P24593
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.