

Product Datasheet

Human IGFBP7 ELISA Kit (orb407420)

Description	Human Insulin-like growth factor-binding protein 7 (IGFBP7/MAC25/PSF) ELISA kit
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Target	IGFBP7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, homogenates
Uniprot ID	Q16270
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.