

## Product Datasheet

### Human AMY2A ELISA Kit (orb407424)

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Description</b>       | Human Pancreatic Amylase,PAMY<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                      |
| <b>Range</b>             | 47 mU/mL-3000 mU/mL                        |
| <b>Target</b>            | AMY2A                                      |
| <b>Note</b>              | For research use only                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm               |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates       |
| <b>Uniprot ID</b>        | <b>P04746</b>                              |
| <b>Sensitivity</b>       | 11.75 mU/mL                                |
| <b>Expiration Date</b>   | Please enquire.                            |