

Product Datasheet

Human INHBB ELISA Kit (orb407453)

Description	human inhibin beta-B,Inhbb ELISA Kit
Reactivity	Human
Range	2 pg/mL-400 pg/mL
Tested Applications	ELISA
Target	INHBB
Note	For research use only
Application notes	Sample Type: serum, plasma Assay Time: 1-5 hours Sample Volume: 50-100 µl Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	P09529
Sensitivity	1 pg/mL
Expiration Date	Please enquire.