

Product Datasheet

Human Cryopyrin ELISA Kit (orb407457)

Description	Human NACHT, LRR and PYD domains-containing protein 3 (NLRP3/C1orf7/CIAS1/NALP3/PYPAF1) ELISA kit
Reactivity	Human
Range	0.156 ng/ml - 10 ng/ml
Target	NLRP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q96P20
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.