

## Product Datasheet

### Human REG4 ELISA Kit (orb407494)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Description</b>       | Human regenerating gene IV,REG-4<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                         |
| <b>Range</b>             | 0.312 ng/mL-20 ng/mL                          |
| <b>Target</b>            | REG4                                          |
| <b>Note</b>              | For research use only                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                  |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates          |
| <b>Uniprot ID</b>        | <b>Q9BYZ8</b>                                 |
| <b>Sensitivity</b>       | 0.078 ng/mL                                   |
| <b>Expiration Date</b>   | Please enquire.                               |