

Product Datasheet

Human REG1A ELISA Kit (orb407495)

Description	Human Regenerating Islet-derived 1 Alpha/pancreatic Stone Protein/pancreatic Thread Protein (REG1A) ...
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	REG1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	P05451
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.