

## Product Datasheet

### Human MUC7 ELISA Kit (orb407505)

|                          |                                                                |
|--------------------------|----------------------------------------------------------------|
| <b>Description</b>       | Human Mucin-7,MUC7 ELISA Kit                                   |
| <b>Reactivity</b>        | Human                                                          |
| <b>Range</b>             | 18.75 ng/mL-1200 ng/mL                                         |
| <b>Target</b>            | MUC7                                                           |
| <b>Note</b>              | For research use only                                          |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                   |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>cell lysates,saliva |
| <b>Uniprot ID</b>        | <b>Q8TAX7</b>                                                  |
| <b>Sensitivity</b>       | 4.68 ng/mL                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                |