

## Product Datasheet

### Human ITGAM ELISA Kit (orb407521)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Human Integrin $\alpha$ M, ITG/CD11b<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                             |
| <b>Range</b>             | 0.78 ng/mL-50 ng/mL                               |
| <b>Target</b>            | ITGAM                                             |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates              |
| <b>Uniprot ID</b>        | <b>P11215</b>                                     |
| <b>Sensitivity</b>       | 0.195 ng/mL                                       |
| <b>Expiration Date</b>   | Please enquire.                                   |