

## Product Datasheet

### Human PNPLA2 ELISA Kit (orb407530)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| <b>Description</b>       | Human Adipose Triglyceride Lipase,ATGL ELISA kit |
| <b>Reactivity</b>        | Human                                            |
| <b>Range</b>             | 6.25 mIU/mL-400 mIU/mL                           |
| <b>Target</b>            | PNPLA2                                           |
| <b>Note</b>              | For research use only                            |
| <b>Application notes</b> | Detection Wavelength: 450 nm                     |
| <b>Sample Types</b>      | serum, urine, tissue homogenates                 |
| <b>Uniprot ID</b>        | <b>Q96AD5</b>                                    |
| <b>Sensitivity</b>       | 1.56 mIU/mL                                      |
| <b>Expiration Date</b>   | Please enquire.                                  |