

## Product Datasheet

### Human LXA4 ELISA Kit (orb407533)

|                          |                                                                                    |
|--------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Lipoxin A4,LXA4 ELISA Kit                                                    |
| <b>Reactivity</b>        | Human                                                                              |
| <b>Range</b>             | 0.625 pg/mL-40 pg/mL                                                               |
| <b>Target</b>            | LXA4                                                                               |
| <b>Note</b>              | For research use only                                                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                       |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>urine, tissue homogenates, cell lysates |
| <b>Sensitivity</b>       | 0.156 pg/mL                                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                                    |