

Product Datasheet

Human TNF ELISA Kit (orb407542)

Description	Human Tumor necrosis factor α , TNF- α ELISA KIT
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	TNF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P01375
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.