

Product Datasheet

Human TNFSF15 ELISA Kit (orb407546)

Description	Human tumor necrosis factor superfamily 15 (TL1A/TNFSF15) ELISA Kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	TNFSF15
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O95150
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.