

Product Datasheet

Human TNFSF9 ELISA Kit (orb407547)

Description	Human Tumor necrosis factor ligand superfamily member 9(TNFSF9) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Target	TNFSF9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P41273
Sensitivity	4.7 pg/mL
Expiration Date	Please enquire.