

Product Datasheet

Human TNFRSF18 ELISA Kit (orb407548)

Description	Human Tumor necrosis factor receptor superfamily member 18(TNFRSF18) ELISA kit
Reactivity	Human
Range	0.156 ng/ml- 10ng/ml
Target	TNFRSF18
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Y5U5
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.