

Product Datasheet

Human TNFRSF6B ELISA Kit (orb407550)

Description	Human Tumor necrosis factor receptor superfamily member 6B(TNFRSF6B) ELISA kit
Reactivity	Human
Range	0.625 ng/ml - 40 ng/ml
Target	TNFRSF6B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O95407
Sensitivity	0.156 ng/ml
Expiration Date	Please enquire.