

Product Datasheet

Human TRAF6 ELISA Kit (orb407552)

Description	Human tumor necrosis factor receptor-associated factor 6, TRAF6 ELISA kit
Reactivity	Human
Range	39 pg/mL-2500 pg/mL
Target	TRAF6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y4K3
Sensitivity	9.8 pg/mL
Expiration Date	Please enquire.