

Product Datasheet

Human TRAIL-R3 ELISA Kit (orb407555)

Description	Human tumor necrosis factor-related apoptosis-inducing ligand receptor 3, TRAIL-R3 ELISA Kit
Reactivity	Human
Range	78 pg/mL-5000 pg/mL
Target	TRAIL-R3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O14798
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.