

Product Datasheet

Human TNFAIP6 ELISA Kit (orb407557)

Description	Human Tumor necrosis factor-inducible gene 6 protein(TNFAIP6) ELISA kit
Reactivity	Human
Range	47 ng/ml - 3000 ng/ml
Target	TNFAIP6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P98066
Sensitivity	11.75 ng/ml
Expiration Date	Please enquire.