

## Product Datasheet

### Human CTSL1 ELISA Kit (orb407590)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Human Cathepsin L1 (CTSL1/CTSL)<br>ELISA kit |
| <b>Reactivity</b>        | Human                                        |
| <b>Range</b>             | 0.312 ng/mL-20 ng/mL                         |
| <b>Target</b>            | CTSL1                                        |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                 |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates         |
| <b>Uniprot ID</b>        | <b>P07711</b>                                |
| <b>Sensitivity</b>       | 0.078 ng/mL                                  |
| <b>Expiration Date</b>   | Please enquire.                              |