

Product Datasheet

Human borrelia burgdorferi antibody (IgM) ELISA Kit (orb407603)

Description	Human borrelia burgdorferi (Bb) antibody (IgM) ELISA kit
Reactivity	Human
Range	Request Information
Target	borrelia burgdorferi (Bb) antibody (IgM)
Note	For research use only
Application notes	Detection Wavelength: 450nm
Sample Types	serum, plasma, cerebrospinal fluid (CSF)
Sensitivity	Request Information
Expiration Date	Please enquire.