

Product Datasheet

Human CXCL17 ELISA Kit (orb407617)

Description	Human chemokine (C-X-C motif) ligand 17 (CXCL17) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	CXCL17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, tissue homogenates
Uniprot ID	Q6UXB2
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.