

Product Datasheet

Human ANGPTL8 ELISA Kit (orb407621)

Description	Human angiopoietin-like protein 8(ANGPTL8) ELISA Kit
Reactivity	Human
Range	6.25 pg/mL-400 pg/mL
Target	ANGPTL8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6UXH0
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.