

Product Datasheet

Monkey Apolipoprotein E ELISA Kit (orb407679)

Description	Monkey Apolipoprotein E(APOE) ELISA kit
Reactivity	Monkey
Range	6.25 ng/mL-400 ng/mL
Target	APOE
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, saliva, tissue homogenates
Uniprot ID	Q28502
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.