

Product Datasheet

Monkey FSH ELISA Kit (orb407690)

Description	Monkey follicle-stimulating hormone,FSH ELISA Kit
Reactivity	Monkey
Range	0.2 mIU/mL-80 mIU/mL
Target	FSH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q3YC02
Sensitivity	0.2 mIU/mL
Expiration Date	Please enquire.