

Product Datasheet

Monkey MCP1 ELISA Kit (orb407691)

Description	Monkey monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP1/MCAF ELISA kit
Reactivity	Monkey
Range	62.5 pg/mL-4000 pg/mL
Target	CCL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P61275
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.