

Product Datasheet

Mouse ADAMTS4 ELISA Kit (orb407717)

Description	Mouse A disintegrin and metalloproteinase with thrombospondin motifs 4(ADAMTS4) ELISA kit
Reactivity	Mouse
Range	3.12 ng/mL-200 ng/ mL
Target	ADAMTS4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8BNJ2
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.