

Product Datasheet

Mouse ADAMTS5 ELISA Kit (orb407718)

Description	Mouse A disintegrin and metalloproteinase with thrombospondin motifs 5(ADAMTS5) ELISA kit
Reactivity	Mouse
Range	0.625 ng/mL-40 ng/mL
Target	ADAMTS5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9R001
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.