

## Product Datasheet

### Mouse ADAMTSL4 ELISA Kit (orb407733)

<b>Description</b>	Mouse ADAMTS-like protein 4(ADAMTSL4) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.56 ng/mL-100 ng/mL
<b>Target</b>	ADAMTSL4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q80T21</a>
<b>Sensitivity</b>	0.39 ng/mL
<b>Expiration Date</b>	Please enquire.