

## Product Datasheet

### Mouse AGRN ELISA Kit (orb407743)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Mouse Agrin(AGRN) ELISA kit          |
| <b>Reactivity</b>        | Mouse                                |
| <b>Range</b>             | 0.312 ng/mL-20 ng/mL                 |
| <b>Target</b>            | AGRN                                 |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450 nm         |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates |
| <b>Uniprot ID</b>        | <b>A2ASQ1</b>                        |
| <b>Sensitivity</b>       | 0.078 ng/mL                          |
| <b>Expiration Date</b>   | Please enquire.                      |